

GENIESORI®

will give you
great happiness
in daily life



(16227)
Gwanggyosan-ro, Yeongtong-gu,
Suwon-si, Gyeonggi-do, Korea
Tel_ +82-31-265-4427
Fax_ +82-50-4373-4428



Personal Sound Amplifier (HA1000)


GENIE invite you to a new world of sound




GENIESORII
Equalizer in My Ear



Causes of hearing loss



Expensive hearing aid price



According to research, listening to loud sounds regardless of age can seriously damage your ears. And maintaining this bad hearing increases the chance of serious accidents in daily life. Now, protect yourself with competitive-cost, high-performance device PSA .

We deliver the best sound quality
for you



UWICOM's Story

"YOUR WISH IS MY COMMAND"

This is the line that Genie, the fairy of Aladdin Lamp, told Aladdin that he can do everything what he want. We have been thinking and studying how to live more comfortable for our family. Someone's wish to hear and communicate better is the reason why we act.



About PSA



01 Automatic fitting algorithm technology








This function uses up to 26 filter channels to create an equalizer that would have the best available audio quality (clean sound without any noise) and inter-band rejection.



02 Best Sound Quality Database

We provide 7 general types of calibrated sound for user. You can select tones easily by just one touch.

* The colors at the bottom are LED lights that show which mode you are listening to.



03 BLUETOOTH SUPPORT

Connectable to all Bluetooth devices
(Ex) Smartphone, TV, SVS-T broadcast)



04 KEY FUNCTIONS

<p>Use audio DSP : best sound quality</p>	<p>Auto-fitting algorithm</p>	<p>Audible frequency band support</p>
<p>Bluetooth™ Version 4.2</p>	<p>Support microphone and Bluetooth sound sources</p>	<p>Lithium rechargeable battery</p>

Use it when you're doing this!



Use in daily life
: When you talk with your family and friends



Using Bluetooth devices (TV, Smartphone, etc)
: When you want to hear outside sounds clearly



Used by a caregiver / doctor
: Easy to communicate with patient




Use of individual agency counselors
: Use to help your staff work smoothly

GENIESORI quality recognized by Korean Public Institutions



- Public Service Centers and Office in Seoul
- National Dementia Center of the entire district
- Ulsan Northern Police Station
- Cheongju District Court
- Incheon Assistive Technology Center
- Branches of Kookmin and Shinhan Bank
- Seoul National University Hospital, Kyung Hee University Medical Center

Product SPEC



Model	GENIESORI2 (HA1000)
Dimensions	28.2(W) X 82(H) X 12.2(T) mm
Weight	27g
Supported Frequency Range	20~20,000Hz
Equalizer Frequency Channel	Max 26 ch
Bluetooth version	4.2
Sound source type	Microphone / Bluetooth input / SVS-T
Continuous operating time	Approximately 9 hours
Battery	Lithium ion polymer 3.7V, 450mAh
Country of origin	Made in Korea
Certifications	KC / CE / FCC / TELEC / ISO9001